

## MRPL48 antibody

### Product Information

Catalog No.:	FNab05336
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	mitochondrial ribosomal protein L48
Synonyms:	Large ribosomal subunit protein mL48 39S ribosomal protein L48, mitochondrial (L48mt, MRP-L48) MRPL48
Observed MW:	20-21 kDa
Uniprot ID :	Q96GC5

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas cancer using FNab05336(MRPL48 antibody) at dilution of 1:50



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab05336(MRPL48 antibody) at dilution of 1:800