

## MRPS2 antibody

### Product Information

Catalog No.:	FNab05351
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	mitochondrial ribosomal protein S2
Synonyms:	Small ribosomal subunit protein uS2m 28S ribosomal protein S2, mitochondrial (MRP-S2, S2mt) MRPS2
Observed MW:	33 kDa
Uniprot ID :	Q9Y399

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human stomach using FNab05351(MRPS2 antibody) at dilution of 1:100



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab05351(MRPS2 antibody) at dilution of 1:3000