

## STK26 antibody

### Product Information


Catalog No.:	FNab05385
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information


Immunogen:	serine/threonine protein kinase MST4
Synonyms:	Serine/threonine-protein kinase 26 MST3 and SOK1-related kinase Mammalian STE20-like protein kinase 4 (MST-4, STE20-like kinase MST4) Serine/threonine-protein kinase MASK STK26 MASK MST4
Observed MW:	47-50 kDa
Uniprot ID :	Q9P289

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human colon cancer using FNab05385(STK26 antibody) at dilution of 1:100



MCF7 cells were subjected to SDS PAGE followed by western blot with FNab05385(STK26 antibody) at dilution of 1:800