

## PWWP3B antibody

### Product Information

Catalog No.:	FNab05438
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	melanoma associated antigen(mutated) 1-like 1
Synonyms:	PWWP domain-containing DNA repair factor 3B (PWWP3B) Mutated melanoma-associated antigen 1-like protein 1 (MUM1-like protein 1) PWWP domain-containing protein MUM1L1 PWWP3B MUM1L1
Observed MW:	None
Uniprot ID :	Q5H9M0

### Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab05438(MUM1L1 antibody) at dilution of 1:50