

## MYLK4 antibody

### Product Information

Catalog No.:	FNab05496
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	myosin light chain kinase family, member 4
Synonyms:	Myosin light chain kinase family member 4 Sugen kinase 85 (SgK085) MYLK4 SGK085
Observed MW:	45 kDa
Uniprot ID :	Q86YV6

### Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart tissue slide using FNab05496(MYLK4 Antibody) at dilution of 1:200



L02 cells were subjected to SDS PAGE followed by western blot with FNab05496(MYLK4 Antibody) at dilution of 1:1000