

NECAB3 antibody

Product Information

Catalog No.:	FNab05640
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibits the interaction of APBA2 with beta-amyloid precursor protein(APP), and hence allows formation of beta-amyloid.

Immunogen information


Immunogen:	N-terminal EF-hand calcium binding protein 3
Synonyms:	N-terminal EF-hand calcium-binding protein 3 Amyloid-beta A4 protein-binding family A member 2-binding protein Nek2-interacting protein 1 Neuronal calcium-binding protein 3 X11L-binding protein 51 NECAB3 APBA2BP NIP1 SYTIP2 XB51
Observed MW:	40-44 kDa
Uniprot ID :	Q96P71

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human pancreas using FNab05640(NECAB3 antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab05640(NECAB3 antibody) at dilution of 1:500