

NRBP2 antibody

Product Information

Catalog No.:	FNab05851
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May regulate apoptosis of neural progenitor cells during their differentiation.

Immunogen information

Immunogen:	nuclear receptor binding protein 2
Synonyms:	Nuclear receptor-binding protein 2 Transformation-related gene 16 protein (TRG-16) NRBP2
Observed MW:	33 kDa, 58 kDa
Uniprot ID :	Q9NSY0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse kidney tissue were subjected to SDS PAGE followed by western blot with FNab05851(NRBP2 antibody) at dilution of 1:500