

NXF3 antibody

Product Information

Catalog No.:	FNab05943
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May function as a tissue-specific nuclear mRNA export factor.

Immunogen information

Immunogen:	nuclear RNA export factor 3
Synonyms:	Nuclear RNA export factor 3 TAP-like protein 3 (TAPL-3) NXF3 TAPL3
Observed MW:	None
Uniprot ID :	Q9H4D5

Application

Reactivity:	Human
Tested Application:	ELISA, IHC
Recommended dilution:	IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas tissue slide using FNab05943(NXF3 antibody) at dilution of 1:50