

OSGEPL1 antibody

Product Information

Catalog No.:	FNab06026
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)


Background

Immunogen information


Immunogen:	O-sialoglycoprotein endopeptidase-like 1
Synonyms:	tRNA N6-adenosine threonylcarbamoyltransferase, mitochondrial N6-L-threonylcarbamoyladenine synthase (t(6)A synthase) O-sialoglycoprotein endopeptidase-like protein 1 (OSGEP-like protein 1) t(6)A37 threonylcarbamoyladenine biosynthesis protein OSGEPL1 tRNA threonylcarbamoyladenine biosynthesis protein OSGEPL1 OSGEPL1 GCP1
Observed MW:	45 kDa
Uniprot ID :	Q9H4B0

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:10 - 1:100
Image:	



Immunohistochemistry of paraffin-embedded human liver injury using FNab06026(OSGEPL1 Antibody) at dilution of 1:100



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06026(OSGEPL1 Antibody) at dilution of 1:1000