

OSMR antibody

Product Information

Catalog No.:	FNab06028
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

Associates with IL31RA to form the IL31 receptor. Binds IL31 to activate STAT3 and possibly STAT1 and STAT5. Capable of transducing OSM-specific signaling events.

Immunogen information

Immunogen:	oncostatin M receptor
Synonyms:	Oncostatin-M-specific receptor subunit beta Interleukin-31 receptor subunit beta (IL-31 receptor subunit beta, IL-31R subunit beta, IL-31R-beta, IL-31RB) OSMR OSMRB
Observed MW:	None
Uniprot ID :	Q99650

Application

Reactivity:	Human, Monkey
Tested Application:	ELISA