

CDK5R1 antibody

Product Information

Catalog No.:	FNab06075
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	cyclin-dependent kinase 5, regulatory subunit 1(p35)
Synonyms:	Cyclin-dependent kinase 5 activator 1 (CDK5 activator 1) Cyclin-dependent kinase 5 regulatory subunit 1 TPKII regulatory subunit Cyclin-dependent kinase 5 activator 1, p35 (p35) Cyclin-dependent kinase 5 activator 1, p25 (p25) Alternative names: Tau protein kinase II 23 kDa subunit (p23) CDK5R1 CDK5R NCK5A
Observed MW:	46 kDa, 35 kDa
Uniprot ID :	Q15078

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:1000; IP: 1:500-1:1000
Image:	



Untreated and Staurosporine treated Jurkat cells were subjected to SDS PAGE followed by western blot with FNab06075(p35 antibody) at dilution of 1:5000