

PLA2G7 antibody

Product Information

Catalog No.:	FNab06107
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	phospholipase A2, group VII(platelet-activating factor acetylhydrolase, plasma)
Synonyms:	Platelet-activating factor acetylhydrolase (PAF acetylhydrolase) 1-alkyl-2-acetyl-glycerophosphocholine esterase 2-acetyl-1-alkyl-glycerophosphocholine esterase Group-VIIA phospholipase A2 (gVIIA-PLA2) LDL-associated phospholipase A2 (LDL-PLA(2)) PAF 2-acylhydrolase PLA2G7 PAFAH
Observed MW:	65 kDa
Uniprot ID :	Q13093

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:200-1:1000; IHC: 1:50-1:200
Image:	



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using FNab06107(PLA2G7 Antibody) at dilution of 1:50



U937 cells were subjected to SDS PAGE followed by western blot with FNab06107(PLA2G7 antibody) at dilution of 1:1000