

PANK1 antibody

Product Information

Catalog No.:	FNab06133
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in the physiological regulation of the intracellular CoA concentration.

Immunogen information

Immunogen:	pantothenate kinase 1
Synonyms:	Pantothenate kinase 1 (hPanK, hPanK1) Pantothenic acid kinase 1 PANK1 PANK
Observed MW:	61-64 kDa, 40-45 kDa
Uniprot ID :	Q8TE04

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



Raji cells were subjected to SDS PAGE followed by western blot with FNab06133(PANK1 antibody) at dilution of 1:300