

PAX3 antibody

Product Information

Catalog No.:	FNab06169
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transcription factor that may regulate cell proliferation, migration and apoptosis. Involved in neural development and myogenesis.

Immunogen information


Immunogen:	paired box 3
Synonyms:	Paired box protein Pax-3 HuP2 PAX3 HUP2
Observed MW:	60 kDa
Uniprot ID :	P23760

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:500-1:1000
Image:	



IP Result of anti-PAX3 (IP: FNab06169, 4ug;
Detection: FNab06169 1:500) with mouse skin
tissue lysate 3000ug.



HEK-293 cells were subjected to SDS PAGE
followed by western blot with FNab06169(PAX3
antibody) at dilution of 1:1000