

## PHLDA1 antibody

### Product Information

Catalog No.:	FNab06396
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Seems to be involved in regulation of apoptosis. May be involved in detachment-mediated programmed cell death. May mediate apoptosis during neuronal development. May be involved in regulation of anti-apoptotic effects of IGF1. May be involved in translational regulation.

### Immunogen information

Immunogen:	pleckstrin homology-like domain, family A, member 1
Synonyms:	Pleckstrin homology-like domain family A member 1 Apoptosis-associated nuclear protein Proline- and glutamine-rich protein (PQ-rich protein, PQR protein) Proline- and histidine-rich protein T-cell death-associated gene 51 protein PHLDA1 PHRIP TDAG51
Observed MW:	40-45 kDa
Uniprot ID :	Q8WV24

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab06396(PHLDA1 antibody) at dilution of 1:500