

PGF antibody

Product Information

Catalog No.:	FNab06440
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Involved in GPI-anchor biosynthesis through the transfer of ethanolamine phosphate to the third mannose of GPI.

Immunogen information


Immunogen:	placental growth factor
Synonyms:	Placenta growth factor (PlGF) PGF PGFL PLGF
Observed MW:	50 kDa
Uniprot ID :	P49763

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:1000; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab06440(PGF Antibody) at dilution of 1:50



human placenta tissue were subjected to SDS PAGE followed by western blot with FNab06440(PGF antibody) at dilution of 1:1000