

EGFLAM antibody

Product Information

Catalog No.:	FNab06454
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	5A5
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	EGF-like, fibronectin type III and laminin G domains
Synonyms:	Pikachurin Agrin-like protein EGF-like, fibronectin type-III and laminin G-like domain-containing protein EGFLAM AGRINL AGRNL PIKA
Observed MW:	111 kDa
Uniprot ID :	Q63HQ2

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



human retina tissue were subjected to SDS PAGE followed by western blot with FNab06454(Pikachurin antibody) at dilution of 1:1000