

PIKFYVE antibody

Product Information

Catalog No.:	FNab06465
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	phosphatidylinositol-3-phosphate/phosphatidylinositol 5-kinase, type III
Synonyms:	1-phosphatidylinositol 3-phosphate 5-kinase (Phosphatidylinositol 3-phosphate 5-kinase) FYVE finger-containing phosphoinositide kinase PIKfyve Phosphatidylinositol 3-phosphate 5-kinase type III (PIPkin-III, Type III PIP kinase) Serine-protein kinase PIKFYVE PIKFYVE KIAA0981 PIP5K3
Observed MW:	237 kDa
Uniprot ID :	Q9Y2I7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse lung tissue were subjected to SDS PAGE followed by western blot with FNab06465(PIP5K3 antibody) at dilution of 1:300