

ALPP antibody

Product Information

Catalog No.:	FNab06520
Size:	100µg
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	2D6
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	alkaline phosphatase, placental(Regan isozyme)
Synonyms:	Alkaline phosphatase, placental type Alkaline phosphatase Regan isozyme Placental alkaline phosphatase 1 (PLAP-1) ALPP PLAP
Observed MW:	66 kDa
Uniprot ID :	P05187

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human placenta tissue slide using FNab06520(ALPP Antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab06520(ALPP Antibody) at dilution of 1:1000