

POU3F3-Specific antibody

Product Information

Catalog No.:	FNab06651
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	POU class 3 homeobox 3
Synonyms:	POU domain, class 3, transcription factor 3 Brain-specific homeobox/POU domain protein 1 (Brain-1, Brn-1) Octamer-binding protein 8 (Oct-8) Octamer-binding transcription factor 8 (OTF-8) POU3F3 BRN1 OTF8
Observed MW:	66 kDa
Uniprot ID :	P20264

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:5000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human brain using FNab06651(POU3F3-Specific antibody) at dilution of 1:100



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab06651(POU3F3-Specific antibody) at dilution of 1:1000