

PRAM1 antibody

Product Information

Catalog No.:	FNab06745
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


May be involved in myeloid differentiation. May be involved in integrin signaling in neutrophils. Binds to PtdIns(4)P.

Immunogen information

Immunogen:	PML-RARA regulated adaptor molecule 1
Synonyms:	PML-RARA-regulated adapter molecule 1 (PRAM, PRAM-1) PRAM1
Observed MW:	70 kDa
Uniprot ID :	Q96QH2

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



U-937 cells were subjected to SDS PAGE followed by western blot with FNab06745(PRAM1 Antibody) at dilution of 1:300