

BBC3 antibody

Product Information

Catalog No.:	FNab06956
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	BCL2 binding component 3
Synonyms:	Bcl-2-binding component 3, isoforms 1/2 JFY-1 p53 up-regulated modulator of apoptosis BBC3 PUMA
Observed MW:	18-21 kDa
Uniprot ID :	Q9BXH1

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IP
Recommended dilution:	WB: 1:500-1:2000; IP: 1:200-1:2000
Image:	



IP Result of anti-PUMA (IP:FNab06956, 4ug;
Detection:FNab06956 1:500) with mouse heart
tissue lysate 3200ug.



A549 cells were subjected to SDS PAGE followed
by western blot with FNab06956(PUMA antibody)
at dilution of 1:500