

## PEX2 antibody

### Product Information


Catalog No.:	FNab06966
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	peroxisomal membrane protein 3, 35kDa
Synonyms:	Peroxisome biogenesis factor 2 35 kDa peroxisomal membrane protein Peroxin-2 Peroxisomal membrane protein 3 Peroxisome assembly factor 1 (PAF-1) RING finger protein 72 PEX2 PAF1 PMP3 PMP35 PXMP3 RNFB72
Observed MW:	35-40 kDa
Uniprot ID :	P28328

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



mouse liver tissue were subjected to SDS PAGE followed by western blot with FNab06966( PXMP3 Antibody) at dilution of 1:600