

RAMP1 antibody

Product Information

Catalog No.:	FNab07097
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Transports the calcitonin gene-related peptide type 1 receptor(CALCRL) to the plasma membrane. Acts as a receptor for calcitonin-gene-related peptide(CGRP) together with CALCRL.

Immunogen information


Immunogen:	receptor(G protein-coupled) activity modifying protein 1
Synonyms:	Receptor activity-modifying protein 1 Calcitonin-receptor-like receptor activity-modifying protein 1 (CRLR activity-modifying protein 1) RAMP1
Observed MW:	17 kDa
Uniprot ID :	O60894

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded rat heart tissue slide using FNab07097(RAMP1 Antibody) at dilution of 1:100



HepG2 cells were subjected to SDS PAGE followed by western blot with FNab07097(RAMP1 antibody) at dilution of 1:1000