

RCAN2 antibody

Product Information

Catalog No.:	FNab07202
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development.

Immunogen information


Immunogen:	regulator of calcineurin 2
Synonyms:	Calcipressin-2 Down syndrome candidate region 1-like 1 Myocyte-enriched calcineurin-interacting protein 2 (MCIP2) Regulator of calcineurin 2 Thyroid hormone-responsive protein ZAKI-4 RCAN2 DSCR1L1 ZAKI4
Observed MW:	24 kDa, 35 kDa
Uniprot ID :	Q14206

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab07202(RCAN2 antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07202(RCAN2 antibody) at dilution of 1:300