

## RELТ antibody

### Product Information

Catalog No.:	FNab07236
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Mediates activation of NF-kappa-B. May play a role in T-cell activation.

### Immunogen information

Immunogen:	RELТ tumor necrosis factor receptor
Synonyms:	Tumor necrosis factor receptor superfamily member 19L Receptor expressed in lymphoid tissues RELТ TNFRSF19L
Observed MW:	47 kDa
Uniprot ID :	Q969Z4

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07236(RELT antibody) at dilution of 1:1000