

## REXO4 antibody

### Product Information

Catalog No.:	FNab07247
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	REX4, RNA exonuclease 4 homolog
Synonyms:	RNA exonuclease 4 Exonuclease XPMC2 Prevents mitotic catastrophe 2 protein homolog (hPMC2) REXO4 PMC2 XPMC2H
Observed MW:	47-50 kDa
Uniprot ID :	Q9GZR2

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human testis using FNab07247(REXO4, REX4 antibody) at dilution of 1:100



Jurkat cells were subjected to SDS PAGE followed by western blot with FNab07247(REXO4, REX4 antibody) at dilution of 1:1000