

RFC2 antibody

Product Information

Catalog No.:	FNab07248
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


The elongation of primed DNA templates by DNA polymerase delta and epsilon requires the action of the accessory proteins proliferating cell nuclear antigen(PCNA) and activator 1. This subunit binds ATP(By similarity).

Immunogen information

Immunogen:	replication factor C(activator 1) 2, 40kDa
Synonyms:	Replication factor C subunit 2 Activator 1 40 kDa subunit (A1 40 kDa subunit) Activator 1 subunit 2 Replication factor C 40 kDa subunit (RF-C 40 kDa subunit, RFC40) RFC2
Observed MW:	35-40 kDa
Uniprot ID :	P35250

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07248(RFC2 Antibody) at dilution of 1:600