

RIMS3 antibody

Product Information

| | |
|---------------|--|
| Catalog No.: | FNab07304 |
| Size: | 100µg |
| Form: | liquid |
| Purification: | Immunogen affinity purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Rabbit |
| Clonality: | polyclonal |
| Clone ID: | None |
| IsoType: | IgG |
| Storage: | PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background


Regulates synaptic membrane exocytosis.

Immunogen information

| | |
|--------------|--|
| Immunogen: | regulating synaptic membrane exocytosis 3 |
| Synonyms: | Regulating synaptic membrane exocytosis protein 3 (Nim3) RIM3 gamma Rab-3-interacting molecule 3 (RIM 3) RIMS3 KIAA0237 |
| Observed MW: | 40 kDa |
| Uniprot ID : | Q9UJD0 |

Application

| | |
|-----------------------|-------------------|
| Reactivity: | Human, Mouse, Rat |
| Tested Application: | ELISA, WB |
| Recommended dilution: | WB: 1:500-1:2000 |
| Image: | |



HL-60 cells were subjected to SDS PAGE followed by western blot with FNab07304(RIMS3 antibody) at dilution of 1:500