

RNASE11 antibody

Product Information

Catalog No.:	FNab07325
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	ribonuclease, RNase A family, 11(non-active)
Synonyms:	Probable ribonuclease 11 (RNase 11) RNASE11 C14orf6
Observed MW:	30 kDa
Uniprot ID :	Q8TAA1

Application

Reactivity:	Human
Tested Application:	ELISA, IHC, WB
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human liver using FNab07325(RNASE11 antibody) at dilution of 1:50



human testis tissue were subjected to SDS PAGE followed by western blot with FNab07325(RNASE11 antibody) at dilution of 1:1000