

## RNPC3 antibody

### Product Information

Catalog No.:	FNab07369
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Participates in pre-mRNA U12-dependent splicing, performed by the minor spliceosome which removes U12-type introns. U12-type introns comprises less than 1% of all non-coding sequences. Binds to the 3'-stem-loop of m(7)G-capped U12 snRNA.

### Immunogen information

Immunogen:	RNA-binding region(RNP1, RRM) containing 3
Synonyms:	RNA-binding region-containing protein 3 RNA-binding motif protein 40 (RNA-binding protein 40) U11/U12 small nuclear ribonucleoprotein 65 kDa protein (U11/U12 snRNP 65 kDa protein, U11/U12-65K) RNPC3 KIAA1839 RBM40 RNP SNRNP65
Observed MW:	62-65 kDa
Uniprot ID :	Q96LT9

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



3T3-L1 cells were subjected to SDS PAGE followed by western blot with FNab07369(RNPC3 antibody) at dilution of 1:600