

RRP1 antibody

Product Information

Catalog No.:	FNab07498
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a critical role in the generation of 28S rRNA.

Immunogen information

Immunogen:	ribosomal RNA processing 1 homolog
Synonyms:	Ribosomal RNA processing protein 1 homolog A Novel nuclear protein 1 (NNP-1) Nucleolar protein Nop52 RRP1-like protein RRP1 D21S2056E NNP1 NOP52 RRP1A
Observed MW:	53 kDa
Uniprot ID :	P56182

Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07498(RRP1 antibody) at dilution of 1:500