

RSRC1 antibody

Product Information

Catalog No.:	FNab07513
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in pre-mRNA splicing. Involved in both constitutive and alternative pre-mRNA splicing. May have a role in the recognition of the 3' splice site during the second step of splicing.

Immunogen information


Immunogen:	arginine/serine-rich coiled-coil 1
Synonyms:	Serine/Arginine-related protein 53 (SRrp53) Arginine/serine-rich coiled-coil protein 1 RSRC1 SRRP53
Observed MW:	38 kDa, 50 kDa
Uniprot ID :	Q96IZ7

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500 - 1:2000; IF: 1:50 - 1:200; IHC: 1:50 - 1:200
Image:	



Immunohistochemistry of paraffin-embedded human rat testis tissue slide using FNab07513(RSRC1 Antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07513(RSRC1 Antibody) at dilution of 1:1000