

SARS2 antibody

Product Information

Catalog No.:	FNab07610
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Catalyzes the attachment of serine to tRNA(Ser). Is also able to aminoacylate tRNA(Sec) with serine, to form the misacylated tRNA L-seryl-tRNA(Sec), which will be further converted into selenocysteinyl-tRNA(Sec)(By similarity).

Immunogen information

Immunogen:	seryl-tRNA synthetase 2, mitochondrial
Synonyms:	Serine--tRNA ligase, mitochondrial SerRSmt Seryl-tRNA synthetase (SerRS) Seryl-tRNA(Ser/Sec) synthetase SARS2 SARSM
Observed MW:	53 kDa
Uniprot ID :	Q9NP81

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500 - 1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07610(SARS2 antibody) at dilution of 1:500