

SEC23IP antibody

Product Information

Catalog No.:	FNab07680
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in the organization of endoplasmic reticulum exit sites. Specifically binds to phosphatidylinositol 3-phosphate(PI(3)P), phosphatidylinositol 4-phosphate(PI(4)P) and phosphatidylinositol 5-phosphate(PI(5)P).

Immunogen information


Immunogen:	SEC23 interacting protein
Synonyms:	SEC23-interacting protein p125 SEC23IP
Observed MW:	125 kDa
Uniprot ID :	Q9Y6Y8

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human heart using FNab07680(SEC23IP antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07680(SEC23IP antibody) at dilution of 1:1000