

EXOC3 antibody

Product Information

Catalog No.:	FNab07686
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information


Immunogen:	exocyst complex component 3
Synonyms:	Exocyst complex component 3 Exocyst complex component Sec6 EXOC3 SEC6 SEC6L1
Observed MW:	86 kDa
Uniprot ID :	O60645

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2400; IF: 1:20-1:200; IHC: 1:50-1:500
Image:	



Immunohistochemistry of paraffin-embedded human gliomas using FNab07686(SEC6 antibody) at dilution of 1:100



human brain tissue were subjected to SDS PAGE followed by western blot with FNab07686(SEC6 antibody) at dilution of 1:86