

anti- Securin / PTTG1 antibody

Product Information

Catalog No.:	FNab07702
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	pituitary tumor-transforming 1
Synonyms:	EAP1, Esp1 associated protein, HPTTG, pituitary tumor transforming 1, PTTG, PTTG1, Securin, Tumor transforming protein 1, TUTR1
Observed MW:	22 kDa
Uniprot ID :	O95997

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, IF
Recommended dilution:IF:	1:10-1:100
Image:	


Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com



Immunofluorescence analysis of HepG2 cells using PTTG1 antibody (FNab07702). Blue: DAPI for nuclear staining.