

SEMA3C antibody

Product Information

Catalog No.:	FNab07706
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Binds to plexin family members and plays an important role in the regulation of developmental processes. Required for normal cardiovascular development during embryogenesis. Functions as attractant for growing axons, and thereby plays an important role in axon growth and axon guidance(By similarity).

Immunogen information

Immunogen:	sema domain, immunoglobulin domain(Ig), short basic domain, secreted,(semaphorin) 3C
Synonyms:	Semaphorin-3C Semaphorin-E (Sema E) SEMA3C SEMAE
Observed MW:	70 kDa
Uniprot ID :	Q99985

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse heart tissue were subjected to SDS PAGE followed by western blot with FNab07706(SEMA3C antibody) at dilution of 1:500