

## SETBP1 antibody

### Product Information

Catalog No.:	FNab07765
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	SET binding protein 1
Synonyms:	SET-binding protein (SEB) SETBP1 KIAA0437
Observed MW:	170-180 kDa
Uniprot ID :	Q9Y6X0

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



HeLa cells were subjected to SDS PAGE followed by western blot with FNab07768( SETBP1 Antibody) at dilution of 1:300