

SFI1 antibody

Product Information

Catalog No.:	FNab07783
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in the dynamic structure of centrosome-associated contractile fibers via its interaction with CETN2.

Immunogen information

Immunogen:	Sfi1 homolog, spindle assembly associated(yeast)
Synonyms:	Protein SFI1 homolog (hSFI1) SFI1 KIAA0542
Observed MW:	148 kDa
Uniprot ID :	A8K8P3

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab07783(SFI1 antibody) at dilution of 1:400