

## SUGP2 antibody

### Product Information

Catalog No.:	FNab07788
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	splicing factor, arginine/serine-rich 14
Synonyms:	SURP and G-patch domain-containing protein 2 Arginine/serine-rich-splicing factor 14 Splicing factor, arginine/serine-rich 14 SUGP2 KIAA0365 SFRS14
Observed MW:	120-130 kDa
Uniprot ID :	Q8IX01

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human prostate cancer using FNab07788(SFRS14 antibody) at dilution of 1:50



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab07788(SFRS14 antibody) at dilution of 1:500