

SHC4 antibody

Product Information

Catalog No.:	FNab07842
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Activates both Ras-dependent and Ras-independent migratory pathways in melanomas.
Contributes to the early phases of agrin-induced tyrosine phosphorylation of CHRNB1.

Immunogen information

Immunogen:	SHC(Src homology 2 domain containing) family, member 4
Synonyms:	SHC-transforming protein 4 Rai-like protein (RaLP) SHC-transforming protein D (hShcD) Src homology 2 domain-containing-transforming protein C4 (SH2 domain protein C4) SHC4 SHCD
Observed MW:	69 kDa
Uniprot ID :	Q6S5L8

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with FNab07842(SHC4 antibody) at dilution of 1:1500