

## SHD antibody

### Product Information

Catalog No.:	FNab07844
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


May function as an adapter protein.

### Immunogen information

Immunogen:	Src homology 2 domain containing transforming protein D
Synonyms:	SH2 domain-containing adapter protein D SHD
Observed MW:	40-50 kDa
Uniprot ID :	Q96IW2

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse brain tissue were subjected to SDS PAGE followed by western blot with FNab07844( SHD Antibody) at dilution of 1:600