

## SKAP1 antibody

### Product Information

Catalog No.:	FNab07891
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information


Immunogen:	src kinase associated phosphoprotein 1
Synonyms:	Src kinase-associated phosphoprotein 1 Src family-associated phosphoprotein 1 Src kinase-associated phosphoprotein of 55 kDa (SKAP-55, pp55) SKAP1 SCAP1 SKAP55
Observed MW:	55 kDa
Uniprot ID :	Q86WV1

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB, IHC
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human spleen tissue slide using FNab07891(SKAP55,SKAP1 Antibody) at dilution of 1:50



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab07891(SKAP55,SKAP1 antibody) at dilution of 1:800