

SLC13A4 antibody

Product Information

Catalog No.:	FNab07909
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Sodium/sulfate cotransporter that mediates sulfate reabsorption in the high endothelial venules(HEV).

Immunogen information

Immunogen:	solute carrier family 13(sodium/sulfate symporters), member 4
Synonyms:	Solute carrier family 13 member 4 Na(+)/sulfate cotransporter SUT-1 NaS2 SLC13A4 SUT1
Observed MW:	126 kDa
Uniprot ID :	Q9UKG4

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



HEK-293 cells were subjected to SDS PAGE followed by western blot with FNab07909(SLC13A4 antibody) at dilution of 1:500