

SLC35F2 antibody

Product Information

Catalog No.:	FNab07954
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Putative solute transporter.

Immunogen information

Immunogen:	solute carrier family 35, member F2
Synonyms:	Solute carrier family 35 member F2 SLC35F2
Observed MW:	41-50 kDa
Uniprot ID :	Q8IXU6

Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



A549 cells were subjected to SDS PAGE followed by western blot with FNab07954(SLC35F2 Antibody) at dilution of 1:600