

GINS4 antibody

Product Information

Catalog No.:	FNab07979
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


None

Immunogen information

Immunogen:	GINS complex subunit 4(Sld5 homolog)
Synonyms:	DNA replication complex GINS protein SLD5 GINS complex subunit 4 DNA replication complex GINS protein SLD5, N-terminally processed GINS4 SLD5
Observed MW:	26 kDa
Uniprot ID :	Q9BRT9

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse thymus tissue were subjected to SDS PAGE followed by western blot with FNab07979(SLD5 antibody) at dilution of 1:300