

SOAT2 antibody

Product Information

Catalog No.:	FNab08098
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background


Plays a role in lipoprotein assembly and dietary cholesterol absorption. In addition to its acyltransferase activity, it may act as a ligase. May provide cholesteryl esters for lipoprotein secretion from hepatocytes and intestinal mucosa.

Immunogen information


Immunogen:	sterol O-acyltransferase 2
Synonyms:	Sterol O-acyltransferase 2 Acyl-coenzyme A:cholesterol acyltransferase 2 (ACAT-2) Cholesterol acyltransferase 2 SOAT2 ACACT2 ACAT2
Observed MW:	46 kDa
Uniprot ID :	O75908

Application

Reactivity:	Human
Tested Application:	ELISA, WB, IHC, IF
Recommended dilution:	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Image:	



Immunohistochemistry of paraffin-embedded human small intestine using FNab08098(SOAT2 antibody) at dilution of 1:50



HeLa cells were subjected to SDS PAGE followed by western blot with FNab08098(SOAT2 antibody) at dilution of 1:1000