

## SOX12 antibody

### Product Information

Catalog No.:	FNab08121
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


Binds to the sequence 5'-AACAAAT-3'.

### Immunogen information

Immunogen:	SRY(sex determining region Y)-box 12
Synonyms:	Transcription factor SOX-12 Protein SOX-22 SOX12 SOX22
Observed MW:	35-40 kDa
Uniprot ID :	O15370

### Application

Reactivity:	Human, Mouse
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



mouse testis tissue were subjected to SDS PAGE followed by western blot with FNab08121( SOX12 Antibody) at dilution of 1:1000