

## SPANXC antibody

### Product Information

Catalog No.:	FNab08151
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Background


None

### Immunogen information

Immunogen:	SPANX family, member C
Synonyms:	Sperm protein associated with the nucleus on the X chromosome C Cancer/testis antigen 11.3 (CT11.3) Cancer/testis-associated protein CTp11 Nuclear-associated protein SPAN-Xc (SPANX-C) SPANX family member C SPANXC
Observed MW:	20 kDa
Uniprot ID :	Q9NY87

### Application

Reactivity:	Human
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:500-1:2000
Image:	



human testis tissue were subjected to SDS PAGE followed by western blot with FNab08151(SPANXC antibody) at dilution of 1:300